# La gamma autoclavi

## Autoclave range

CE  
2460

### Caratteristiche tecniche Technical characteristics

**ONYX B 5.0**  
carichi tipo B  
type B loads

**ONYX B 7.0**  
carichi tipo B  
type B loads

Dimensioni (LxAxP) / Dimensions (WxHxD)	475 x 495 x 650 mm	474 X 497 X 720 mm
Peso / Weight	60 Kg	64 Kg
Tensione/ Power supply	230 V - 50/60 HZ	230 V - 50/60 HZ
Assorbimento/ Absorption	1,5 kW	2,2 kW

### Dimensione camera • Chamber dimension

Diametro camera / Chamber diameter	Ø 245 mm	Ø 245 mm
Profondità / Depth	320 mm	460 mm

### Optional

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EN

Autoclave is the key device within each medical structure, therefore it is fundamental choosing a well-known and referenced product which can rely on a widespread network of technical support. Onyx B is the new range of Tecno-Gaz autoclaves, a world leading Company for many years.

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## Double load... Same time

Using standard size autoclave about 4 kg instruments can be sterilized, while with Onyx B you can sterilize from 5 kg to 7 kg of instruments. **With 1 double load you are saving 50% of your time.**

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### SD card

SD CARD system allows **constant machine updating** by a simple download on it and then by inserting it into the autoclave system. **An extraordinary option in order to have constantly efficient and updated device.**

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### Cycles for all needs

Flexibility is fundamental to rationalise working and create operative efficiency. **Onyx B autoclaves offer a wide range of cycles all studied to help you every day.**

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### Best management of instruments

**Having big load autoclave helps you to rationalize instruments flow and increases your operative efficiency.** Efficiency means profit, Onyx B autoclaves are advantageous for this too.

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### Energy saving

Onyx B autoclaves have low energetic consumption as well as unbeatable consumption/ load ratio. **Economic saving and benefit for the environment.**

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EN

### Big load

New regulations and new sterilization concept have developed sterilization practice in all medical and dental facilities along with wider range of instruments and spread of professional facilities which developed the need of autoclaves able to sterilize large loads. **Onyx project is born with the aim to satisfy these need adopting specific and new solutions enabling to manage easily also partial loads.**

### FL System - Full Load System

To avoid manual loads with possibility of water pouring off, Onyx B autoclaves throughout their automatic clean water tank automatic system are using a level sensor control pump with automatic lock. **Frontal loading, no opening cover or removing system, enabling simple placement of the machine.**

### Proper trays for any type of load

Trays of different lengths enables you to insert instruments inside bulky envelopes, to place load more rationally enabling sterilization of **different loads following the correct guidelines for load placement.**

### PHD System

Preheating system with stand-by for minimum temperature **maintenance and consequent overall sterilization cycles time reduction.**

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### Water net connection

The autoclaves Onyx B are autonomous as equipped with double tank of pure water and exhausted water, they are prepared for the connection to the water supply through osmosis system mod. PURA Plus.

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### Inspectionable tanks

In case of insertion of unsuitable liquids, autoclave tanks are inspectionable and cleanable as per technical normative regulations. Onyx B autoclaves are supplied with pure and used water tanks. Each tank is supplied with probes indicating levels.

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### Overnight cycles

Automatic overnight cycles with anti-condensation system and automatic switch off.

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### ID system

Authentication intergrated system of the operator which is protected by password and permits at the end of the cycle to associate operator to sterilized load having visual convalidation of the operator himself. The same authentication is also reported by the procedure of load traceability.

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### Reduced maintenance

Onyx B autoclaves allow you to make big loads thus reducing the number of cycles. This means extending autoclave life over 50% reducing maintenance costs and amortizing the cost of the product over the time.

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EN

## 7 sterilisation cycles

Onyx B is born with the aim of rationalizing the work of sterilization operators, therefore the autoclaves are equipped with cycles which have been expressively conceived for improving time management and operator's needs within outpatient clinics.

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
## Completely automatic cycles

All the cycles are completely automatic, **all the phases are managed by the advanced card equipped with powerful microprocessor**. Each and every phase (vacuum phase, preheating, exposure time, drying) is monitored, backed up and managed automatically. There is no human involvement in the carrying out of the cycles.

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EN

Europa B Evo is an autoclave that has been developed jointly by the Tecno-Gaz research office and important external development facilities. Europa B Evo is indispensable for all medical and dental facilities and in all fields where safe sterilization of all materials is required: hollow, porous, solid, wrapped and unwrapped. Europa B Evo continues the success that Tecno-Gaz has achieved over the years, in fact, thousands of professionals use Tecno-Gaz autoclaves and are very satisfied.

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### P.I.D. System

P.I.D. means Proportional Integral Derivate and represents the apex of closed loop controls. The structure of this technologically advanced control is simple and flexible. P.I.D. ensures a constant control and monitoring of the whole operating system of the autoclave, acting automatically in case of anomalies or incorrect values in technical and electronic parameters. **This device behaves like a "real technician" with real-time diagnostic and intervention capability. This solution reduces anomalies and warnings, granting maximum reliability and minimising the machine inactivity times.**

### Safety on the door

We have applied three safety devices on the autoclave door:

- **The first device is a mechanic and electromagnetic one:** it triggers when the sterilization cycle begins and completely blocks the door opening. This device can be removed at the end of the cycle by pushing the device-unlocking button, but only if positive pressure is near zero.
- **The second device is an electrical one:** it triggers automatically when an attempt to force the door open is detected.
- **The third device is dynamic** and acts on the structure of the gasket.

### Safety devices in place

- Certified safety valve.
- Ped on sterilization chambers.
- Insulation of the chamber with certified ceramic material.

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### Barometric alignment

**Automatic barometric alignment** is provided from the first time the device is switched on, without any technical intervention.

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### Overnight cycles

**Automatic overnight cycles** with anti-condensing system and automatic switch-off, with memory effect.

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### Provision for mains water supply

**The osmosis system mod. Pura is an extremely compact** piece of equipment for osmosis, created to be directly connected to the autoclave, guaranteeing a constant supply.

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### Built-in printer

**Built-in printer for the registration of cycles carried out. The printer includes:**

- date, time
- cycle number
- serial number of the device
- minute by minute cycle data
- outcome of the cycle

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EN

The Pura Plus osmosis system is extremely compact and it's designed for connection directly to the autoclave, to guarantee a constant supply of pure water. The system can be installed in a number of positions to facilitate mounting and the optimum organisation of the available space. Pura Plus is fitted with a series of osmotic membranes and filters, with an average duration of about 800 cycles, depending on the original hardness of the water. The cartridge should therefore need replacing about once a year, or after 800 cycles. The overall cost of the filters is very reasonable, and the cost per litre is therefore almost negligible. When the Pura Plus system is not producing water, it is not pressurized, and this avoids all kinds of problems caused by the accidental breakage of couplings and fitting.

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## Why an autoclave osmosis system?

The osmosis system is recognised as the best method, ensuring treated water free from impurities and hard particles in order to protect the workings of the autoclave. The system allows for a direct connection with the autoclave, ensuring service and a reduction in costs relating to the purchase of distilled water.

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**SIR Oftalmica Srl**