



STERIBEST™ Trays


- Designed to protect delicate instruments during sterilization
- Silicone pin mats included for added protection
- Made from the polymer Ultem® to ensure strength and long lasting durability
- Transparent lid for easy identification of contents
- Trays are stackable

ENHANCED FEATURES:


- **Raised node base surface** and **slotted lid** facilitate circulation of steam and aid in drying
- Extra slots mean increased steam penetration and more rapid drying
- Helps eliminate wet pack




Item **CP614** *Base, Lid, Mat*
6 in. x 2.5 in. x .75 in.
152mm x 64mm x 19mm




Item **CP778** *Base, Lid, Mat*
7.5 in. x 2.5 in. x .75 in.
190.5mm x 63.5mm x 19mm




Item **CP471D1** *Base, Lid, Mat*
7 in. x 4 in. x 1.5 in.
178mm x 102mm x 38mm
(not pictured)



Item **CP471S** *Base, Lid, Mat*
7 in. x 4 in. x .75 in.
178mm x 102mm x 19mm
(not pictured)



Item **CP634** *Base, Lid, Mat*
6.5 in. x 4.0 in. x .75 in.
165mm x 102mm x 19mm



Item **CP1038D2**
Base, Lid, Mat, Insert Tray, Insert Mat
10 in. x 6.0 in. x 1.5 in.
254mm x 152mm x 38mm

Item **CP1038D1** *Base, Lid, Mat*
10 in. x 6.0 in. x 1.5 in.
254mm x 152mm x 38mm
(not pictured)

Item **CP1038S** *Base, Lid, Mat*
10 in. x 6.0 in. x .75 in.
254mm x 152mm x 19mm
(not pictured)



Manufacturer of SteriBest Trays
RICA SURGICAL PRODUCTS, INC. Since 1984
9207 Ivanhoe Street | Schiller Park, IL 60176
800-889-3218 | 847-671-5600 | Fax 847-671-5690
www.ricasurgical.com

* Note: All dimensions are I.D.
Instruments shown for illustrative purposes
and are not included with the instrument trays.