




Suggested Anwar Keratoplasty Set

1-218

Anwar Corneal Scissor - Right




- Blunt tips, curved blades
- Long inner blade with smoothly polished ends
- Cut length 5.0mm, blade length 11.0mm
- Round handle, length 108.0mm



1-219

Anwar Corneal Scissor - Left




- Blunt tips, curved blades
- Long outer blade with smoothly polished ends
- Cut length 5.0mm, blade length 11.0mm
- Round handle, length 108.0mm



6-112

Anwar Keratoplasty Hook




- 0.18mm diameter tip for 1.0mm length
- Flat face at tip
- 45° angled shaft, tip to angle length 10.0mm
- Round handle, length 118.0mm



6-099-3



Anwar Keratoplasty Spatula



- 0.25 diameter blunt tip tapered from 0.5mm diameter
- Angled shaft, tip to angle length 7.0mm
- Round handle, length 116.0mm




9-735-1
Anwar Keratoplasty 16 Blade Ring Marker



- 16 radial marks
- 10.5mm outside diameter ring
- Peripheral cut-outs spaced around the blades
- Central crossbar
- Round handle, length 125mm

DK7003
Sterilising Tray



- External dimensions: 17.5cm x 11.0cm x 2.6cm
- Internal dimensions: 16.5cm x 10.0cm x 1.8cm (without silicone mat)
- One silicone mat

www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.com



© March 2011 Duckworth & Kent

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

at the Leading Edge