



SMART CP[®] 11
Auto Chart Projector



Simple Operation and Various Visual Chart



Wireless Control

By operating MEDIZS Application, wireless communication allows users to control both Chart Projector and Phoropter. Reduced step of controlling visual chart provides fast measurement to both users and patients.



Easy Tightening Knob

The easy tightening knob allows a simple installation to the unit table stand. (No screws) / Patent pending



Easy Maintenance

Convenient assembly of bulb screw cap helps easy maintenance and replacement. / Patent No.KR10-1273450



Power Indication LED

An elegant white LED backlight indicates status of the device.



Selection of Lamp

Lamp can be selected by the user's preference, either replaceable LED or halogen, depending on the measurement environment.



One-touch Open

To open the case of the instrument, users simply press the two buttons simultaneously instead of using any tool.

CP-11 Specifications

Chart / 41 Selections

Mask	3 horizontal lines 5 vertical lines 21 single letters
Chart Selection Speed	Average 0.2 sec
Filter	Red/Green, Polarizing filter
Projection Distance	2.0m - 7.0m (3m is recommendable)
Projection Magnification	30X (at 5m)
Lamp	Halogen(12v 30W) / P-LED(3W)
Power Saving	3min / 5min / 10min

Type

Smart CP-11	RS-232 Communication
Smart CP-11B	Bluetooth Communication
Chart Type	A type (Fraction indication) B type (Feet indication) C type (Meter indication) D type (For Russian)

Power

Power Supply	AC 100-240V, 50/60Hz (Free voltage)
Power Consumption	Halogen: 85-105VA / LED: 35-45VA

Size

Dimension (mm)	Main Body: 204[W] X 313[D] X 260[H] Remote controller: 64[W] X 195[D] X 20[H]
Net Weight	Main body_ 3.8kg / 8.4lbs Remote controller_ 180g / 0.41lbs

Above specification and design may be altered without notice.

MEDIZS Inc.	TEL, +82,42,933,5531 / FAX, +82,42,933,5525 94-10, Techno 2-ro, Yuseong-gu, Daejeon, S,Korea, 305-509 Email : medizs@medizs.com
MEDIZS USA, Inc.	TEL, +1,310-212-3444 / FAX, +1,310-212-3446 2441 W, 205th Street Suite C204 Torrance, CA 90503 Email : support@medizs.com