

# SMART RT<sup>®</sup> 13

## Refraction Table Series



RT-13R



RT-13S



## RT-13 Refraction Unit Table

Made with a compact design for a stylish optometry room, it is the smartest unit table that can change the height of a chair using iPad.



### NEW Controller

Operator can reduce the refraction time by easily controlling the device using App via Bluetooth wireless network.



### Ergonomic Design

Motorized vertical adjustment up to 350mm provide convenience for measuring position regardless of patient's body size.



### Clean Environment

- Horizontal Sliding Table helps convenient exam even in a small room.
- Low noise level provides clean refraction environment when handling table.



### Storage space

Wide space helps to store trial lens set.

### RT-13 Specifications

Power Supply	AC 220V 50/60Hz
Lamp	25W light bulb
Table Rotating Angle	90°
Arm Angle	34°
Chair Move	350mm(Max)
Size	Total_ 1070(W)*1125(D)*1731(H)mm Table_ 1070(W)*490(D)*30(T)mm Moving Table_ 285(W)*375(D)*20(T)mm
Weight	Unit_ 170kg Chair_ 60kg



## RT-13R Compact Refraction Table



### Ergonomic Design

- Up and Down function helps measurement regardless of patient's body size.
- Low noise level provides clean refraction environment when handling table.



### Arm Rotation

Arm rotation angle until 90 degree helps to reduce inconvenience significantly before and after exam.



### Storage space

Wide space helps to store trial lens set.

### RT-13R Specifications

Power Supply	AC 220V 60Hz
Lamp	25W light bulb
Arm Rotating Angle	90°
Size	Table_700(W)*650(D)mm Up&Down_660~925mm Arm_640(H)mm
Weight	57kg

## RT-13S Single Table

### RT-13S Specifications

Power Supply	AC 220V 50/60Hz
Power Consumption	40VA
Size	Table_560(W)*490(D)mm Up&Down_695~995mm
Weight	22kg (Allowable weight: 70kg)

### MEDIZS Inc.

TEL, +82,42,933,5531 / FAX, +82,42,933,5525  
94-10, Techno 2-ro, Yuseong-gu, Daejeon, S.Korea, 305-509  
Email : medizs@medizs.com

### MEDIZS USA Inc.

TEL : +1,310,212,3444 / FAX : +1,310,212,3446  
2441 W. 205th Street Suite C204 Torrance, CA 90503  
Email : support@medizs.com

### DISTRIBUTED by.