

# VeraPlug™ FlexFit™

Punctal Occlusion System

The VeraPlug™ FlexFit™ is designed for **simple sizing, easy insertion, patient comfort, and retention**—providing a proven treatment for chronic dry eye.

- Two sizes fit most patients for **simple sizing**
- FlexFit™ nose design provides **easy insertion** and **superior retention**
- Low-profile dome for excellent **patient comfort**

### 30-DAY RETENTION PLEDGE\*

If a VeraPlug™ FlexFit™ is lost within 30 days of insertion, Lacriversa will replace it free of charge.



#### AVAILABLE CONFIGURATIONS

##### Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

| SIZE    | 1 PAIR   |
|---------|----------|
| X-Small | VFF-5000 |
| Small   | VFF-5001 |
| Medium  | VFF-5002 |
| Large   | VFF-5003 |

##### Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 10 pairs per box. Inserters included.

| SIZE    | 10 PAIRS |
|---------|----------|
| X-Small | BFF-5000 |
| Small   | BFF-5001 |
| Medium  | BFF-5002 |
| Large   | BFF-5003 |



#### VERAPLUG™ FLEXFIT™ SIZE CHART

| PUNCTAL OPENING | SIZE    |
|-----------------|---------|
| 0.2mm to 0.3mm  | X-Small |
| 0.3mm to 0.5mm  | Small   |
| 0.6mm to 0.8mm  | Medium  |
| 0.9mm to 1.0mm  | Large   |

# VeraPlug™

## Punctal Occlusion System

The **VeraPlug™** is designed to provide a simple and effective treatment for chronic dry eye.

- Engineered to perfectly fit the puncta for **superior retention**
- Low-profile dome for **enhanced patient comfort**
- Proprietary shaft design for **easy insertion** and **proper anatomic fit**
- Available in four sizes



Choose the **sterile preloaded** configuration for the ultimate in convenience or **nonsterile bulk** for increased value.

### SIZE CHART

| PUNCTAL OPENING | VERAPLUG™ SIZE |
|-----------------|----------------|
| 0.4mm to 0.6mm  | Small          |
| 0.6mm to 0.7mm  | Medium         |
| 0.7mm to 0.8mm  | Large          |
| 0.8mm to 1.0mm  | X-Large        |



### AVAILABLE CONFIGURATIONS

#### Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

| SIZE    | 1 PAIR   |
|---------|----------|
| Small   | VSP-1001 |
| Medium  | VSP-1002 |
| Large   | VSP-1003 |
| X-Large | VSP-1004 |

#### Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 10 pairs per box. Inserters included.

| SIZE    | 10 PAIRS |
|---------|----------|
| Small   | BNS-1001 |
| Medium  | BNS-1002 |
| Large   | BNS-1003 |
| X-Large | BNS-1004 |

# VeraPlug™ flow

Punctal Occlusion System



The **VeraPlug™ Flow** is designed to provide a simple and effective means of slowing tear drainage.

- Inner channel reduces but does not eliminate tear flow, providing only **partial occlusion**
- Engineered to perfectly fit the puncta for **superior retention**
- Low-profile dome for **enhanced patient comfort**
- Proprietary shaft design for **easy insertion** and **proper anatomic fit**

## AVAILABLE CONFIGURATIONS

### Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

| SIZE    | 1 PAIR   |
|---------|----------|
| Small   | VFL-1001 |
| Medium  | VFL-1002 |
| Large   | VFL-1003 |
| X-Large | VFL-1004 |

### Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 5 pairs per box. Inserters included.

| SIZE    | 5 PAIRS  |
|---------|----------|
| Small   | BVF-1001 |
| Medium  | BVF-1002 |
| Large   | BVF-1003 |
| X-Large | BVF-1004 |

Choose the **sterile preloaded** configuration for the ultimate in convenience or **nonsterile bulk** for increased value.



## SIZE CHART

| PUNCTAL OPENING | VERAPLUG™ FLOW SIZE |
|-----------------|---------------------|
| 0.4mm to 0.6mm  | Small               |
| 0.6mm to 0.7mm  | Medium              |
| 0.7mm to 0.8mm  | Large               |
| 0.8mm to 1.0mm  | X-Large             |

# Vera180™

Synthetic Absorbable Lacrimal Plugs

Vera180™ synthetic absorbable lacrimal plugs provide approximately 6 months of treatment for dry eye syndrome.

- Provides extended temporary occlusion lasting approximately **180 days**
- Ideal for treatment of **post-surgical dry eye**, and to determine potential effectiveness of a non-absorbable punctal plug
- Also indicated for use in treating the dry eye components of various **ocular surface diseases**, including contact lens intolerance




**6**  
Months  
OF RELIEF FROM DRY EYE

## Vera180™ SINGLE PAIR

Synthetic absorbable lacrimal plugs made of poly-p-dioxanone (PDO). Absorbs in approximately 180 days. Two plugs per box.

SIZE

|       |           |
|-------|-----------|
| 0.2mm | ESK-1002S |
| 0.3mm | ESK-1003S |
| 0.4mm | ESK-1004S |
| 0.5mm | ESK-1005S |

## Vera180™ TEN PAIRS

Synthetic absorbable lacrimal plugs made of poly-p-dioxanone (PDO). Absorbs in approximately 180 days. 20 plugs per box (10 packs of 2 plugs each).

SIZE

|       |          |
|-------|----------|
| 0.2mm | ESK-1002 |
| 0.3mm | ESK-1003 |
| 0.4mm | ESK-1004 |
| 0.5mm | ESK-1005 |

# VeraC7™

Collagen Inserts

# Vera90™

Synthetic Extended Inserts



**VeraC7™** and **Vera90™** are intended for temporary punctal occlusion following surgery or as a trial for permanent punctal occlusion in dry eye therapy.

## VeraC7™

Collagen punctal inserts for diagnostic use, providing highly effective occlusion in the first 2 to 3 days post-insertion. Absorbs completely in 7 to 10 days. 60 inserts per box (10 packs of 6 inserts each).

### SIZE

|       |          |
|-------|----------|
| 0.2mm | CSP-1002 |
| 0.3mm | CSP-1003 |
| 0.4mm | CSP-1004 |

- Provides temporary occlusion lasting **7 to 10 days**
- Ideal for **diagnostic use** in chronic dry eye patients considering permanent punctal occlusion
- Also indicated for extending medication contact time and effectiveness



## Vera90™


Synthetic extended punctal inserts made of E-Caprolactone-L-Lactide copolymer (PCL). Absorbs in 60 to 180 days. 20 inserts per box (10 packs of 2 inserts each).

### SIZE

|       |          |
|-------|----------|
| 0.2mm | ESP-1002 |
| 0.3mm | ESP-1003 |
| 0.4mm | ESP-1004 |
| 0.5mm | ESP-1005 |

- Provides extended temporary occlusion lasting **60 to 180 days**
- Ideal for the treatment of **dry eye following LASIK** or other surgical procedures
- Also indicated for use in treating the dry eye components of various ocular surface diseases





**SIR** oftalmica

DISTRIBUTING EXCELLENCE SINCE 1989

SIR OFTALMICA SPA

VIA NAPOLEONA 7

22100 – COMO

ITALY

TEL. +39.031.570869

E-MAIL: [INFO@SIROFTALMICA.COM](mailto:INFO@SIROFTALMICA.COM)

[WWW.SIROFTALMICA.COM](http://WWW.SIROFTALMICA.COM)