

VeraPlug™ FlexFit™

Punctal Occlusion System

The VeraPlug™ FlexFit™ is designed for **simple sizing, easy insertion, patient comfort, and retention**—providing a proven treatment for chronic dry eye.

- Two sizes fit most patients for **simple sizing**
- FlexFit™ nose design provides **easy insertion** and **superior retention**
- Low-profile dome for excellent **patient comfort**

30-DAY RETENTION PLEDGE*

If a VeraPlug™ FlexFit™ is lost within 30 days of insertion, Lacriversa will replace it free of charge.



AVAILABLE CONFIGURATIONS

Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

SIZE	1 PAIR
X-Small	VFF-5000
Small	VFF-5001
Medium	VFF-5002
Large	VFF-5003

Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 10 pairs per box. Inserters included.

SIZE	10 PAIRS
X-Small	BFF-5000
Small	BFF-5001
Medium	BFF-5002
Large	BFF-5003



VERAPLUG™ FLEXFIT™ SIZE CHART

PUNCTAL OPENING	SIZE
0.2mm to 0.3mm	X-Small
0.3mm to 0.5mm	Small
0.6mm to 0.8mm	Medium
0.9mm to 1.0mm	Large

VeraPlug™

Punctal Occlusion System

The **VeraPlug™** is designed to provide a simple and effective treatment for chronic dry eye.

- Engineered to perfectly fit the puncta for **superior retention**
- Low-profile dome for **enhanced patient comfort**
- Proprietary shaft design for **easy insertion** and **proper anatomic fit**
- Available in four sizes



Choose the **sterile preloaded** configuration for the ultimate in convenience or **nonsterile bulk** for increased value.

SIZE CHART

PUNCTAL OPENING	VERAPLUG™ SIZE
0.4mm to 0.6mm	Small
0.6mm to 0.7mm	Medium
0.7mm to 0.8mm	Large
0.8mm to 1.0mm	X-Large



AVAILABLE CONFIGURATIONS

Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

SIZE	1 PAIR
Small	VSP-1001
Medium	VSP-1002
Large	VSP-1003
X-Large	VSP-1004

Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 10 pairs per box. Inserters included.

SIZE	10 PAIRS
Small	BNS-1001
Medium	BNS-1002
Large	BNS-1003
X-Large	BNS-1004

VeraPlug™

Punctal Occlusion System

The **VeraPlug™** is designed to provide a simple and effective treatment for chronic dry eye.

- Engineered to perfectly fit the puncta for **superior retention**
- Low-profile dome for **enhanced patient comfort**
- Proprietary shaft design for **easy insertion** and **proper anatomic fit**
- Available in four sizes



Choose the **sterile preloaded** configuration for the ultimate in convenience or **nonsterile bulk** for increased value.

SIZE CHART

PUNCTAL OPENING	VERAPLUG™ SIZE
0.4mm to 0.6mm	Small
0.6mm to 0.7mm	Medium
0.7mm to 0.8mm	Large
0.8mm to 1.0mm	X-Large



AVAILABLE CONFIGURATIONS

Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

SIZE	1 PAIR
Small	VSP-1001
Medium	VSP-1002
Large	VSP-1003
X-Large	VSP-1004

Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 10 pairs per box. Inserters included.

SIZE	10 PAIRS
Small	BNS-1001
Medium	BNS-1002
Large	BNS-1003
X-Large	BNS-1004

VeraPlug™ flow

Punctal Occlusion System



The **VeraPlug™ Flow** is designed to provide a simple and effective means of slowing tear drainage.

- Inner channel reduces but does not eliminate tear flow, providing only **partial occlusion**
- Engineered to perfectly fit the puncta for **superior retention**
- Low-profile dome for **enhanced patient comfort**
- Proprietary shaft design for **easy insertion** and **proper anatomic fit**

AVAILABLE CONFIGURATIONS

Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

SIZE	1 PAIR
Small	VFL-1001
Medium	VFL-1002
Large	VFL-1003
X-Large	VFL-1004

Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 5 pairs per box. Inserters included.

SIZE	5 PAIRS
Small	BVF-1001
Medium	BVF-1002
Large	BVF-1003
X-Large	BVF-1004

Choose the **sterile preloaded** configuration for the ultimate in convenience or **nonsterile bulk** for increased value.



SIZE CHART

PUNCTAL OPENING	VERAPLUG™ FLOW SIZE
0.4mm to 0.6mm	Small
0.6mm to 0.7mm	Medium
0.7mm to 0.8mm	Large
0.8mm to 1.0mm	X-Large

Vera180™

Synthetic Absorbable Lacrimal Plugs

Vera180™ synthetic absorbable lacrimal plugs provide approximately 6 months of treatment for dry eye syndrome.

- Provides extended temporary occlusion lasting approximately **180 days**
- Ideal for treatment of **post-surgical dry eye**, and to determine potential effectiveness of a non-absorbable punctal plug
- Also indicated for use in treating the dry eye components of various **ocular surface diseases**, including contact lens intolerance



6
Months
OF RELIEF FROM DRY EYE

Vera180™ SINGLE PAIR

Synthetic absorbable lacrimal plugs made of poly-p-dioxanone (PDO). Absorbs in approximately 180 days. Two plugs per box.

SIZE

0.2mm	ESK-1002S
0.3mm	ESK-1003S
0.4mm	ESK-1004S
0.5mm	ESK-1005S

Vera180™ TEN PAIRS

Synthetic absorbable lacrimal plugs made of poly-p-dioxanone (PDO). Absorbs in approximately 180 days. 20 plugs per box (10 packs of 2 plugs each).

SIZE

0.2mm	ESK-1002
0.3mm	ESK-1003
0.4mm	ESK-1004
0.5mm	ESK-1005

VeraC7™

Collagen Inserts

Vera90™

Synthetic Extended Inserts



VeraC7™ and **Vera90™** are intended for temporary punctal occlusion following surgery or as a trial for permanent punctal occlusion in dry eye therapy.

VeraC7™

Collagen punctal inserts for diagnostic use, providing highly effective occlusion in the first 2 to 3 days post-insertion. Absorbs completely in 7 to 10 days. 60 inserts per box (10 packs of 6 inserts each).

SIZE

0.2mm	CSP-1002
0.3mm	CSP-1003
0.4mm	CSP-1004

- Provides temporary occlusion lasting **7 to 10 days**
- Ideal for **diagnostic use** in chronic dry eye patients considering permanent punctal occlusion
- Also indicated for extending medication contact time and effectiveness



Vera90™


Synthetic extended punctal inserts made of E-Caprolactone-L-Lactide copolymer (PCL). Absorbs in 60 to 180 days. 20 inserts per box (10 packs of 2 inserts each).

SIZE

0.2mm	ESP-1002
0.3mm	ESP-1003
0.4mm	ESP-1004
0.5mm	ESP-1005

- Provides extended temporary occlusion lasting **60 to 180 days**
- Ideal for the treatment of **dry eye following LASIK** or other surgical procedures
- Also indicated for use in treating the dry eye components of various ocular surface diseases





SIR oftalmica

DISTRIBUTING EXCELLENCE SINCE 1989

SIR OFTALMICA SPA

VIA NAPOLEONA 7

22100 – COMO

ITALY

TEL. +39.031.570869

E-MAIL: INFO@SIROFTALMICA.COM

WWW.SIROFTALMICA.COM